# MORPHEUS LT

ANAESTHESIA  
WORKSTATION  
FOR ADULTS  
CHILDREN  
AND NEWBORNS

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**Siare**

# MORPHEUS LT

## ANAESTHESIA WORKSTATION WORKSTATION PER ANESTESIA

The Morpheus LT is an anaesthesia workstation suitable for treatment of adults, children and new-borns.

It is provided with a mechanical 3-gas flowmeters group with double scale for low flows and with a group for rapid connection of two vaporizers in line and it is supplied with an advanced manual ventilation system, APL system, or AUX system with traditional To-and-Fro valve switched on and off automatically by ventilator.

The unit frame is made for 30% of light alloy and for 70% of plastic material.

Thanks to the material used, the Morpheus has a total weight of less than 80 kgs.

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The Morpheus LT is equipped with 3 wide and large drawers and a large work-shelf to facilitate operating rooms procedures.

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**SICRE**

## ELECTRONIC VENTILATOR

The ventilator works by a microprocessor electronic control system with a simple and intuitive interface, which allows the anesthetist to operate in a clear and safe way. The system is provided with 3 dedicated softwares to have automatically the default settings for adult and paediatric patients and new-borns.

It is provided with the automatic compensation of dead space, the leakages calculation and the automatic calibration of  $\text{FiO}_2$ .

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PAW  
RATE  
VM - VTE  
T. INSP  
 $\text{FiO}_2$

## PATIENT CIRCUIT GROUP

The patient group of Morpheus LT is compact, lightweight for manual removal and fully autoclavable.

The patient circuit group is supplied with an automatic heating system avoiding humidity with possibility to use IS PAC disposable  $\text{CO}_2$  absorbers of 0.5 kgs or reusable and rechargeable  $\text{CO}_2$  absorbers (Top special of 1 kg).

Both absorbers are with quick-release enabling replacement also during anaesthesia.

The group is supplied with a bellow that can be manually inspected.

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## RESPIRATORY MONITOR MODULE

The respiratory monitor module measures the respiratory parameters that are shown on the TFT 6,5" colour display.

The data reading is made through the reusable flow sensor that incorporates the innovative SIARE technology based on the magnetic oscillation granting an accurate reading of patient respiratory data.

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Pressione; Flusso; Volume; Loop (Volume/Pressione; Volume/Flusso)

## DATA FOR ORDERING

### MORPHEUS LT ANAESTHESIA WORKSTATION

#### OM3.A3s - Morpheus LT

- Anaesthesia Workstation including:
- Mixing device at 3 gas O<sub>2</sub>-N<sub>2</sub>O-Air.
- Electronic ventilator.
- Patient group for ventilation in open, semi-closed, closed at low flow system.
- Support for the rapid connection of two vaporizers.
- TOP SPECIAL CO<sub>2</sub> absorber.

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**Code 100145/IDF**  
Sevoflurane Vaporizer

**Code 100151/IDF**  
Desflurane Vaporizer

**Code 06751**  
Stainless steel patient circuit supporting arm

**Code 713011 CT**  
TFT 6,5" respiratory colour monitor

**Code 30005/10**  
IS PAC  
0,5 Kg CO<sub>2</sub> absorber disposable  
BOX of 10 pcs

**Code 30900**  
TOP SPECIAL  
1 Kg CO<sub>2</sub> absorber rechargeable

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€ 0476

*Made in Italy*

SIARE is not only research, designing and manufacturing but with an efficient and expertise technical service around the world grants the maximum efficiency of its products.

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