EC-3H/ABP COMBINED HOLTER SYSTEM (BLUETOOTH)

- FULL-DISCLOSURE AMBULATORY ECG AND BP RECORDS
- 3 INDEPENDENT ECG CHANNELS
- OSCILLOMETRIC BP MEASUREMENT
- PACEMAKER CHANNEL AND ANALYSIS
- PATIENT EVENT BUTTON
- RECHARGEABLE BATTERIES FOR **ECONOMIC OPERATION**
- USER-FRIENDLY SOFTWARE
- FAST ORS CLASSIFICATION AND RHYTHM ANALYSIS
- HRV, ST, QT, PM, ABP AND AF **ANALYSIS**
- COLOR-CODED GRAPHS AND TABLES
- VARIOUS EDITING AND PRINTING **FUNCTIONS**









EC-3H/ABP

The EC-3H/ABP combines a 3-channel ECG Holter and an ambulatory blood pressure monitor. The ORS and PM can be detected with high precision, while template and rhythm analysis can be interactively modified by the user. This user-friendly software provides high-quality time- and frequency-based HRV, QT, ST, PM, ABP and AF results. The transparency and the interpretation of the results are aided by color-coded graphs and tables. The software offers a wide range of languages to choose from.





TABLE VIEW



FULL VIEW + BP



MAIL:

www:

LABTECH LTD., 4. VÁG STREET DEBRECEN 4031 HUNGARY HTTP://WWW.LABTECH.HU

PHONE: +36-52-310-128 +36-52-412-023 FAX: EMAIL: MEDICAL@LABTECH.HU