

# OXI-CAPNOGRAPHY

Comdekk®

## Portable Oxi-Capnography

MD-660P is a portable multifunctional Capnography monitor providing SpO<sub>2</sub>, pulse, EtCO<sub>2</sub>, Respiration with accuracy and reliable performance.

It features with light and large LCM display providing graphical SpO<sub>2</sub>, EtCO<sub>2</sub> and Trend function on all in one design. Moreover, through the optional software of Smart Link and USB connection to a computer, user can search, assemble and analysis all of patient's information with 72 hours memory. This model uses a side stream sampling system to measure the CO<sub>2</sub> concentration in the patient's breath and respiration frequency and the alarm is preset for unusual oxygen saturation and pulse rate.

MD-660P is an essential part for every emergency & anesthesiology system, and sports center...etc.

