



Elevate the Way You Operate

The CUDA HD Wireless Camera illuminates your procedure more than ever before while also capturing a clear, accurate view for the whole team.

The CUDA LHD-30 Wireless HD Headlight Camera System allows users to capture the surgeon's unique point of view of surgical procedures with unlimited mobility and superior light output. The high-quality camera provides the entire operating room team with a crisp, clear image of the surgical site.

The CUDA LHD-30 can be used for a wide range of applications including improving patient records, refining techniques, documenting revisions and teaching purposes.

Wireless transmission from the recording unit/battery pack to the receiver eliminates additional cables in the operating room and allows the user more comfortable movement around the patient table.

Videos and still images are captured easily using a remote programming device or by pressing the corresponding buttons on the recording unit/battery pack.

Recordings and pictures are saved to a 128G SD card internal to the recording unit and exported easily to a computer via a USB to mini-USB cable.

A robust, one of a kind system designed to provide the surgeon and the OR Team with a unique perspective during surgical procedures without ever compromising mobility, quality of light or light output!



Specifications

Wireless HD Headlight Camera System

- Captures crisp, high-definition video and images
- Control video recording and still images with recording unit or remote programming device
- Real-time image view on the OR Display
- Easy transfer of files to PC
- Optional lenses available for different fields of view
- 70,000 Lux (*measured at 16 in/ 40 cm working distance*) with defined uniform white light throughout the variable spot size
- Battery-powered unit allows maximum mobility
- Lightweight Headband with Rear Cranial Support
- Robust Aluminum Module with variable 20mm to 110mm spot
- System Equipped with Extended Linkage, which accommodates most surgical loupes
- Removable camera for headlight-only use



CUDA[®]
SURGICAL

6018 Bowdendale Avenue, Jacksonville, Florida 32216
Ph: 877-814-2237 • Fx: (+1) 904-733-4832
www.cuda.com

MADE IN THE USA

CE



ISO 13485 Certified